



2025.06.14. <DNA Gel Excision>

- Agarose gel (LE agarose 2.25 g + 1X TAE Buffer 150 ml + SYBR green 4 ul)
- Electrophoresis at 100 V for 20 min.
- PCR product = 16 ul, 1kb plus DNA Ladder = 10 ul



확대샷



1kb plus
DNA Ladder RP RP

13-1 13-2 13-3 14-1 14-2 14-3 15-1 15-2 15-3 16-1 16-2 16-3

RP RP

17-1 17-2 17-3 18-1 18-2 18-3 19-1 19-2 19-3 20-1 20-2 20-3 21-1

RP P.C. P.C. 1kb plus

21-2 21-3 22-1 22-2 22-3 23-1 23-2 23-3 24-1 24-2 24-3 N.C. 1 2 DNA Ladder